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Allen He

Product Manager-consumables & reagents

June 2015

## Milestones & Millipore Firsts

**1954 - Millipore founded (first US based membrane company)**

**1956 - Millipore filter certified by a state laboratory for coliform analysis**

**1964 - First ultrafiltration membrane developed**

**1967 - First centrifugal ultrafiltration (**Centriflo**) product developed**

**1971 - First sterile syringe filter (**Millex®**)**

**1975 - Southern hybridization (**HA membrane**) described**

**1979 - Western blotting, (**HA membrane**) described - Towbin.**

**1979 - **Durapore®** membrane manufactured.**

**1983 - Introduction of first centrifugal ultrafilter with “reverse spin” (**Centricon®**)**

**1985 - Introduction of PVDF blotting membrane (Immobilon™ P)**

**1996 - Introduction of low-binding, high flow membrane (Express)**

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







- ◆ Pore size (孔径)
- ◆ Hydrophilic / Hydrophobic (亲/疏水性)
- ◆ Compatibility (化学相容性)
- ◆ Liquid and Air Flow Rates (流速)
- ◆ Porosity (成孔率)
- ◆ Bubble Point (泡点)
- ◆ Protein Binding (蛋白结合)

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

## F.D.A. Definition:

“A sterilizing filter is one which, when challenged with the Microorganism *B. diminuta*, at a minimum concentration of  $10^7$  organisms per  $\text{cm}^2$  of filter surface, will produce a sterile effluent.”

Pore Size	Common Particles Removed
0.1 $\mu\text{m}$	Mycoplasma, large viruses
0.22 $\mu\text{m}$	USP Sterility certified for bacterial removal ( <i>B. diminuta</i> )
0.45 $\mu\text{m}$	Bacteria ( <i>E.coli</i> , <i>S. aureus</i> , etc)
0.8 $\mu\text{m}$	Molds, yeast
1.0 $\mu\text{m}$	Beads, small mammalian cells
3-5 $\mu\text{m}$	Blood cells, lymphocytes, radioactive particles in air
5-10 $\mu\text{m}$	Mammalian cells
10-40 $\mu\text{m}$	Machined parts contaminants
40-60 $\mu\text{m}$	Tissue, cell aggregates, precipitates
60-180 $\mu\text{m}$	Whole organisms (parasites)

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

## Binding assays 结合实验

- Receptor Ligand Binding (RLB): Glass FC or B
- Kinase: Glass or P-81

## Immunoassays 免疫分析

- Luminex Beads (BV or HV)
- ELISPOT


## Sample Preparation 样本制备

- Neonatal Screening (HV) 新生儿筛查
- Genomics (PCR/ Seq /Dye Terminator Removal)
- Solvinert
- 96 Well Chromatography
- Nylon Mesh


# MultiScreen HTS 96/384-Well Filter Plates

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

# The MultiScreen HTS Manifold = Key Accessory for Multiscreen Plates

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


MultiScreen Well 中的HA or IP膜

# Neonatal Screening(新生儿筛查)

MultiScreen-HV plates are a high throughput method for reproducibly preparing newborn blood samples for PKU(苯丙酮尿症) screening

**A published study showed that MultiScreen plates improve sample preparation efficiency by almost 2x.**

*Screening 3 (1995) 209–223*

- Approximately 1 out of every 16,000 babies is born with phenylketonuria (PKU)
- PKU is an inherited disease caused by the absence of an enzyme for the amino acid phenylalanine
- If not treated, PKU can lead to severe mental retardation and brain damage
- MAHV plates (.45um)

# MultiScreen PCR<sub>μ96</sub> Plates

1. Optimized purification of PCR reaction volumes less than 150 μL
2. A one piece, automation compatible
3. **A format enabling smaller resuspension volumes (20μL)**
4. The ability to reduce PCR reagent usage through its micro well design
5. Higher concentration factors (up to 7.5x) than our standard PCR 96 product
6. Greater than 99% primer removal
7. **Better recovery of smaller fragments (137 bp)**

## PCR<sub>96</sub> 96-well Plate

MSNU03010  
MSNU03050

LSKMPCR10  
LSKMPCR50


# Standard vs. Micro

	<ul style="list-style-type: none"> <li>• <b>PCR<sub>μ96</sub></b></li> <li>• (LSKMPCR10 / 50)</li> </ul>	<ul style="list-style-type: none"> <li>• <b>Standard PCR<sub>96</sub></b></li> <li>• (MSNU03010 / 50)</li> </ul>
<b>Design</b>	<ul style="list-style-type: none"> <li>• One piece design</li> </ul>	<ul style="list-style-type: none"> <li>• 2-piece / underdrain</li> </ul>
<b>Reaction volume input</b>	<ul style="list-style-type: none"> <li>• 1 - 150 μl</li> </ul>	<ul style="list-style-type: none"> <li>• 150 - 300 μl</li> </ul>
<b>Minimum sample recovery volume</b>	<ul style="list-style-type: none"> <li>• 20 μl</li> </ul>	<ul style="list-style-type: none"> <li>• 50 μl</li> </ul>
<b>Automation</b>	<ul style="list-style-type: none"> <li>• Optimal</li> </ul>	<ul style="list-style-type: none"> <li>• Optimal</li> </ul>
<b>Small PCR fragment recovery</b>	<ul style="list-style-type: none"> <li>• &gt;80%</li> </ul>	<ul style="list-style-type: none"> <li>• &gt;50%</li> </ul>

# PAMPA plate – uses the ‘Classic’ design

Part #'s:  
 Standard  
 No underdrain  
 Disp. Receiver  
 Teflon Receiver

MAIPN4550  
 MAIPNTR10  
 MATRNPS50  
 MSSACCEPT0R


# Solvinert & Deepwell Solvinert Plates

MultiScreen 96-well filter plate for organic solvent compatible sample preparation.

## Featuring:

- Choice of 0.5mL or **1.9 mL** versions
- Hydrophilic and hydrophobic 0.45 $\mu$ m high-throughput **PTFE** membranes
- Hydrophobic with optional pre-filter (Deep well only)
- Patented one-piece automation friendly design,
- Broad solvent compatibility and low extractables
- Currently recommended for use with Affymetrix SNP 6.0 prep kit (MDRLN0410)

Drug Discovery and Development often requires the use of organic solvents, strong acids, and other chemicals not compatible with standard filter plates

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
Merck Millipore is a division of MERCK

Multiply your success.  
Multiwell solutions for discovery  
research and sample prep.

Merck Millipore is a division of MERCK

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# What is the Microcon DNA Fast Flow?

Microcon DNA Fast Flow is optimized for the concentration and recovery of gDNA. The low non-specific binding characteristics of the membrane coupled with its unique inverted spin recovery allows the device to provide consistency and reproducibility for gDNA applications.

## Features:

- High recovery for small volumes with reverse spin
- Low binding membrane
- Fast processing time
- Flat membrane design for;
  - Forensic gDNA applications
  - Collection of filtrate

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Purify biologically active proteins with gentle, all-in-one recovery.

Amicon® Pro Purification System

Merck Millipore is a division of Merck KGaA, Darmstadt, Germany

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

# Millicell® Cell Culture Inserts:


## *Advantages of Cells Grown on Filters*

Cells grown on filters experience:

- *Improved cell morphology*
- *Better cell differentiation*
- *More intracellular organelles*
- *Higher cell densities*

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**MF-Millipore™ MCE Membrane  
(mixed cellulose esters)**

For exceptional anatomical and functional polarization.

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**Isopore™ PCF Membrane  
(polycarbonate)**

For growth of attachment-dependent cells without matrix.

**Biopore™ PFTE Membrane  
(polytetrafluoroethylene)**

For low protein binding, live cell viewing, and immunofluorescent applications.

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**PET Membrane  
(polyethylene terephthalate)**

For growth of attachment-dependent cells without matrix.

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
96-well cell culture plate, single-well feeder tray, 96-well receiver tray, and lid	PCF	0.4 μm	1	PSHT004R1
	PET	1 μm		PSRP004R1
96-well cell culture plate, 96-well receiver tray, and lid	PCF	0.4 μm	5	PSHT004S5
96-well cell culture plate, single-well feeder tray, and lid	PCF	0.4 μm	5	PSHT004R5
	PET	1 μm		PSRP004R5

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## Chemotaxis

Directed cell movement toward chemical gradient

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## Haptotaxis

Directed cell movement toward ECM protein

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## Invasion

Cell invasion through an ECM protein &/or another cell layer

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	<b>Millicell</b> 站立式 (孔径)	<b>Millicell</b> 悬挂式 (孔径)	<b>Millicell</b> 24-孔板 (孔径)	<b>Millicell</b> 96-孔板 (孔径)
血管发生	PCF (3, 8)	PET (3, 5, 8)	PET (3, 5, 8)	MultiScreen MIC Plate*
细胞增殖	PCF (0.4)	PET (1)	PCF (0.4)	PCF (0.4) PET (1)
细胞表面受体	PCF (0.4) HA (0.45) CM (0.4)	PET (1)	PCF (0.4) PET (1.0)	PCF (0.4) PET (1)
趋化性	PCF (3, 8)	PET (3, 5, 8)	PET (3, 5, 8)	MultiScreen MIC Plate*
共培养	PCF (0.4) CM (1)	PET (1, 3)	PET (1) PCF (0.4)	PCF (0.4) PET (1)
侵袭性	PCF (8)	PET (5, 8)	PET (5, 8)	MultiScreen MIC Plate*
上皮细胞生长	PET (1)	PET (1)	PET (1)	PCF (0.4) PET (1)
饲养层	PCF (0.4, 3, 8)	PET (0.4, 1, 3, 5, 8)	PET (0.4, 3, 5, 8)	PCF (0.4) PET (1)
荧光检测/免疫组织化学	PCF CM (0.4)	PET (all)	PCF PET (1)	PCF PET (1)
体外受精	CM (0.4)	PET (1)	PET (1)	PET (1)
体外毒理学	PCF (0.4) CM (0.4) HA (0.45)	PET (0.4, 1)	PCF (0.4) PET (1)	PCF (0.4) PET (1)
微生物附属物	PCF (0.4) CM (0.4) HA (0.45)	PET (0.4, 1)	PCF (0.4) PET (1)	PCF (0.4) PET (1)
器官型	Organotypic (0.4)			
相差显微镜检查	CM (0.4)	PET (1)	PET (1)	PET (1)
极化蛋白质分泌物	PCF (0.4) CM (1)	PET (0.4, 1)	PCF (0.4) PET (1)	PCF (0.4) PET (1)
极化摄取	PCF (0.4) CM (0.4) HA (0.45)	PET (0.4, 1)	PCF (0.4) PET (1)	PCF (0.4) PET (1)
转运/渗透性	PCF (0.4)	PET (1)	PCF (0.4) PET (1)	PCF (0.4) PET (1)
肿瘤细胞转移和侵袭	PCF (8)	PET (5, 8)	PET (5, 8)	MultiScreen MIC Plate*

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Description	Plate material/color	Qty/Pk	Sterile	Catalogue No.
MultiScreen® 8-well strip with Immobilon®-P membrane	Acrylic/Clear	10	Yes	M8IPS4510
MultiScreen® <sub>HTS</sub> -IP plate with Immobilon®-P membrane	Acrylic/White	10	Yes	MSIPS4W10
MultiScreen® <sub>HTS</sub> -IP plate with Immobilon®-P membrane	Acrylic/Clear	10	Yes	MSIPS4510
MultiScreen® <sub>HTS</sub> plate with Immobilon®-P membrane without underdrain	Acrylic/White	10	Yes	MAIPSWU10
MultiScreen®-IP plate with Immobilon®-P membrane without underdrain	Acrylic/Clear	10	Yes	MAIPS4510
MultiScreen®-HA with MCE membrane	Styrene/Clear	10	Yes	MAHAS4510

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